



ENDOTRANS MOTORISED CART FOR ENDOSCOPY

FEATURES

- Custom design service provides you with a unique tailored product.
- Robust stylish design to complement the equipment carried.
- Some of the options that can be included on the Endo-Trans are: lockable drawers, endoscope trays, lids, shelves, CO₂ bottle holder, stent holder, eye level organiser with illumination and monitor arm.
- Motorisation available.

Part number

104200

Description

A wide variety of trolleys are available in different formats and specifications. Please contact us to discuss your requirements in more detail.

SCOPE HANGERS

Available in a variety of designs and manufactured from tough durable ABS plastic.

FEATURES

- Hygienic and easy to clean.
- Design provides maximum support for the umbilical.
- Supplied complete with wall brackets, fixing screws and wall plugs.



Part number

103657

Description

Flush wall mounting for GI scopes

Pack size

1

103658

90° mounting for GI scopes

1

103659

Overhead mounting for GI scopes

1

103663

Flush wall mounting for bronchoscopes

1



MONITOR TROLLEY

A specialist trolley suited for a wide range of monitors.

FEATURES

- Will support a wide range of LCD monitors.
- Cable management.
- Ergonomically designed and easy to clean.

Part number

104190

Description

Monitor trolley

Pack size

1



ARTICULATING ARMS

Articulating arms available in a variety of designs to support current range of LCD monitors.

FEATURES

- Easily adjustable.
- Cable management.
- Easily cleaned.
- A wide variety of monitor mounts are available in different styles and specifications. Please contact us to discuss your requirements in more detail.

MEDICAL ISOLATION TRANSFORMER

Designed and built to comply with the requirements of Medical Equipment Specification, EN 60601-1

FEATURES

- Available in a variety of designs. Please call for details.
- Various power options ranging from 750VA to 3000VA, in Toroidal or Laminated construction.
- Full medical use compliance.
- Fused input and output (+e-) providing ultra-safe operation.
- Variety of output couplings can be supplied and enclosures can be custom built to fit customer requirements.

